

## Materials

- Streptavidin coated magnetic beads (e.g., Pierce Cat #PI88816)
- 1.5 mL microcentrifuge tubes
- SiaFind<sup>™</sup> Kit(s), Biotinylated; one or more of the following:
  - Pan-Specific Lectenz<sup>®</sup>: 10 mg/mL stock (Lectenz Bio Cat #SP0502B)
  - $\alpha$ 2,3-Specific Lectenz<sup>®</sup>: 10 mg/mL stock (Lectenz Bio Cat #SP2302B)
  - $\alpha$ 2,6-Specific Reagent: 1 mg/mL stock (Lectenz Bio Cat #SP2602B)
- Binding Buffer:
  - 5X SiaFind<sup>™</sup> Binding Buffer 1 (Lectenz Bio Cat #BA0101) diluted to 1X (SBB1) in Milli-Q water for Pan-Specific and  $\alpha$ 2,3-Specific Lectenz<sup>®</sup>
  - 5X SiaFind<sup>™</sup> Binding Buffer 2 (Lectenz Bio Cat #BA0102) diluted to 1X (SBB2) in Milli-Q water for  $\alpha$ 2,6-Specific Reagent
- Wash Buffer:
  - 5X SiaFind<sup>™</sup> Binding Buffer 1 (Lectenz Bio Cat #BA0101) diluted to 1X + 0.1% Tween-20 (SBB1T)
  - 5X SiaFind<sup>™</sup> Binding Buffer 2 (Lectenz Bio Cat #BA0102) diluted to 1X + 0.1% Tween-20 (SBB2T)
- Elution Buffer:
  - 5X SiaFind<sup>™</sup> Regeneration Buffer 1 (Lectenz Bio Cat #BA0301) diluted to 1X in Milli-Q water (SRB1)
- Alternate elution buffer: SDS-PAGE reducing sample buffer
- Antigen sample (clarified lysates) diluted in SBB1 or SBB2
- Magnetic stand (e.g., DynaMag<sup>™</sup>-2 Magnet, Invitrogen Cat #12321D)

## Pre-washing Streptavidin Magnetic Beads

**Note:** To ensure homogeneity, mix the beads thoroughly before use by repeated inversion or vortexing.

1. Add 0.5 mg (50  $\mu$ L) of streptavidin magnetic beads to a 1.5 mL microcentrifuge tube.
2. Place the tube in the magnetic stand to collect the beads against the side of the tube. Remove and discard the supernatant.
3. Add 1 mL of Wash Buffer to the tube. Invert the tube several times or vortex gently to mix. Collect the beads with a magnetic stand, then remove and discard the supernatant. Repeat 3 times total.

**Note:** Do not allow beads to dry. If necessary, store beads in binding buffer before proceeding with the immunoprecipitation protocol.

## Immunoprecipitation with downstream analysis by non-biotinylated analyte

**Note:** This protocol will allow for greater precipitation yields, but if your antigen sample contains multiple glycosylation sites, some precomplexed biotinylated SiaFind<sup>™</sup> will co-elute. This will show up if subsequent analysis is done using a streptavidin conjugate.

1. Pre-complexing: Combine the 500  $\mu$ L diluted antigen sample (the analyte) with 50  $\mu$ g of biotinylated SiaFind<sup>™</sup> reagent (the “antibody”). Incubate 2 h or overnight at 4°C with mixing.

**Note:** Dilute each sample to a minimum volume of 300-500  $\mu\text{L}$  with cell lysis buffer or Binding/Wash Buffer. We recommend  $\sim 50 \mu\text{g}$  of antigen sample to 50  $\mu\text{g}$  SiaFind™ reagent as a starting point.

2. Add the antigen sample/biotinylated antibody pre-complex to a 1.5 mL microcentrifuge tube containing the pre-washed magnetic beads, followed by incubation for 30 min at room temperature or 1 h at 4°C with mixing.
3. Collect the beads with a magnetic stand and remove and save the supernatant for downstream analysis by SDS-PAGE.
4. Add 500  $\mu\text{L}$  of Wash Buffer to the tube and gently mix. Collect the beads and then collect the supernatant for downstream analysis by SDS-PAGE. Repeat this wash twice.

### Immunoprecipitation with downstream analysis by biotinylated analyte

**Note:** This protocol allows for lower yields, but removes potential SiaFind™ reagent contaminants from the eluate. This is useful if using a streptavidin-conjugate in downstream analysis.

1. Combine the 50  $\mu\text{L}$  streptavidin magnetic beads and 50  $\mu\text{g}$  biotinylated SiaFind™ reagent (the “antibody”) and incubate 30 min at room temperature or 1 h at 4°C with mixing.

**Note:** Dilute each sample to a minimum volume of 300-500  $\mu\text{L}$  with Binding/Wash Buffer.

2. Place the tube into the magnetic stand to collect the beads against the side of the tube. Remove and discard the supernatant.
3. Add 1 mL of Wash Buffer to the tube. Invert the tube several times or vortex gently to mix. Collect the beads with a magnetic stand, then remove and discard the supernatant. Repeat 3 times total.
4. Add the 500  $\mu\text{L}$  antigen sample and incubate 2 h at room temperature or overnight at 4°C with mixing.
5. Collect the beads with a magnetic stand and remove and save the supernatant for analysis.
6. Add 500  $\mu\text{L}$  of Wash Buffer to the tube and gently mix. Collect the beads and then collect the supernatant downstream analysis by SDS-PAGE. Repeat this wash twice.

### Elution Buffer Recovery of Antigen:

1. Add 100  $\mu\text{L}$  of Elution Buffer to the tube. Incubate 10 min at room temperature with mixing.
2. Magnetically separate the beads and save the supernatant containing target antigen.

**Note:** If the eluate is only to be used for SDS-PAGE analysis, the elution can also be done by adding SDS sample buffer to the beads and boiling for 10-15 min at 95°C.

### Downstream analysis:

Downstream analysis typically consists of SDS-PAGE followed by Coomassie blue staining and/or Western blotting against the antigen of interest.