Rev. 2025-06-17

β-N-Acetylglucosaminidase contents

Catalog #	Description	Size	M. W.	Purity	pH range	Storage
GE1201	β- <i>N</i> -Acetylglucosaminidase	500 units, lyophilized	101,594	> 95%	3.5-7.5	-20°C, up to 12 months
BA0801	10X Reaction Buffer 4	1 mL			7.5	4 to 25°C
BA1101	10X Reaction Buffer 7	1 mL			5.0	4 to 25°C

This product is for research use only and not for resale or for any use in the manufacture of a therapeutic or for any diagnostic purpose

Product description: This product is recombinant β -N-Acetylglucosaminidase (glycosyl hydrolase family GH20; E.C. 3.2.1.52), cloned from *Streptococcus pneumoniae* and expressed in *Escherichia coli* with an N-terminal 8xHis tag. The 8xHis tag may be removed by digestion with FasTEVTM (Cat #GE0501), a TEV protease with enhanced stability and catalytic activity.

This enzyme catalyzes the hydrolysis of the non-reducing terminal *N*-Acetylglucosamine (GlcNAc) from oligosaccharides and glycoprotein substrates.



This product does not contain any detectable activities of proteases or other glycosidases.

β-N-Acetylglucosaminidase is supplied with two 10X Reaction Buffers to ensure optimal digestion and ease of use. Reaction Buffer 4 (Cat #BA0801) is used for reactions that require higher than neutral pH buffering, and Reaction Buffer 7 (Cat #BA1101) is the optimal buffer for most digestions.

Unit definition: One unit is defined as the amount of enzyme required to catalyze the release of 1 nmole of p-nitrophenol (pNP) from p-nitrophenyl-N-acetyl- β -D-glucosaminide (pNP-GlcNAc) in 1 min at 37°C in 100 μ L 1X Reaction Buffer 4 (50 mM Tris-HCl, 100 mM NaCl, pH 7.5).

Activity assay: One unit of enzyme is added to $100~\mu L$ of $500~\mu M$ pNP-GlcNAc in 1X Reaction Buffer 7 (50~mM sodium citrate, pH 5.0) at $37^{\circ}C$ for 30~min, followed by addition of $100~\mu L$ of a stop solution (0.2~M sodium borate, pH 9.8). Measure absorption at 405~nm on a plate reader.

Product reconstitution: Dissolve the lyophilized product in 100 μ L molecular grade water to make a 5,000 units/mL (Cat #GE1201) solution in 1X Reaction Buffer 4. Once reconstituted, store at 4°C for up to 7 days or -20°C for up to 3 months. Aliquoting is recommended to avoid repeated freeze-thaw cycles.

Suggested protocol for removal of GlcNAc from glycoproteins:

1. Mix the following components in a microfuge tube:

Glycoprotein (e.g. G0N-BSA; Cat #NG0101) 10 μ g 10X Reaction Buffer 7 (Cat #BA1101) 2 μ L β -N-Acetylglucosaminidase (Cat #GE1201) 1 μ L (5 units) Molecular grade water to 20 μ L final volume

- 2. Incubate at 37°C for 1h.
- 3. Analyze 2 μ L reaction by Western/lectin blot (or ELISA) to determine the completion of reaction. Suggested 1° probe: lectin GSL II, biotinylated.

Reference: Clarke VA, et al. J Biol Chem. 1995 Apr 14;270(15):8805-14. PMID: 7721787